
## Bodycad announces US FDA 510k clearance for BC Fine Osteotomy ™

April 12, 2020

Bodycad was granted FDA 510(k) clearance for the BC Fine Osteotomy<sup>™</sup> system, the first personalized planning and plating system for osteotomies around the knee.

BC Fine Osteotomy™ is synonymous with our highly accurate Bodycad proprietary software, offering a complete solution adapted to each unique patient. It includes a 3D printed patient-specific surgical guide (approved by the surgeon) and a patient specific plate with a small titanium wedge to confirm the 3D angular correction. The BC Fine Osteotomy™ system includes calibrated drill guides with positive stops to protect posterior neurovascular structures, multi-planar 3D planning to increase accuracy and is delivered as a procedure in a box to drive efficiency in the OR.



## Contact

Bodycad USA Inc.

5935 Barclay Lane

Naples, FL 34110

Steve Pruitt VP of US Sales – mobile: 404.307.9974

Stephen Herrington President of US and Asia - mobile: 239.877.7069

info@bodycad.com