

FOR IMMEDIATE RELEASE

CellRight Technologies Announces Launch of Matrix OI® Cortical Fibers

04.05.2017

Jesus Hernandez

President & CEO

Jhernandez@cellrighttechnologies.com

UNIVERSAL CITY, Texas (April 5, 2017) CellRight Technologies, the leader in the development, manufacture and distribution of revolutionary osteoinductive matrices, is pleased to announce the availability of Matrix OI® Cortical Fibers as the latest *in vivo* verified osteoinductive matrix to its growing portfolio of regenerative human derived osteobiologics.

CellRight is recognized in the industry for verifying osteoinductivity post-sterilization in the *in vivo* nude rat model for all product lots. Matrix OI® Cortical Fibers is also tested *in vitro* for the endogenous BMP-2 test marker for osteoinductivity. "CellRight Technologies continues to raise the bar higher by its commitment to provide clinicians and their patients with verified osteoinductive products. Matrix OI Cortical Fibers is available freeze-dried, providing surgeons with an opportunity to rehydrate the Matrix OI Cortical Fibers into a malleable inductive grafting material using the patient's bone marrow aspirate (BMA), blood and/or other cellular concentrates in accordance with the physician's well-informed medical judgment," said Jesus Hernandez, Chief Executive Officer for CellRight Technologies.

Matrix OI® Cortical Fibers are provided sterile (Sterile "R") with a sterility assurance level (SAL) 10^{-6} . Matrix OI® Cortical Fibers have a shelf life of 5 years, do not contain any extrinsic carriers and are made from 100% human bone. Once rehydrated, Matrix OI® Cortical Fibers become malleable with exceptional handling characteristics that resist irrigation. The hydrated Matrix OI® Cortical Fibers may be loaded in a syringe and have the ability to be delivered through an open cannula for use in minimally invasive surgical procedures. The product is available prepackaged in jars for easy mixing in volumes of 1, 2.5, 5 and 10 ccs.

Matrix OI® Cortical Fibers are indicated for use for the treatment of surgically created or traumatic skeletal defects. Matrix OI® Cortical Fibers are regulated as a minimally manipulated product.

For further information, contact CellRight Technologies Customer Sales at (210) 659.9353 or customercare@cellrighttechnologies.com