

Bodycad Gains Regulatory Approval in Canada for BC Fine Osteotomy



May 20, 2020

Bodycad was approved by Health Canada for its BC Fine Osteotomy™ product, the first personalized planning and plating system for osteotomies around the knee.

BC Fine Osteotomy™, similar to Bodycad proprietary software, offers a procedural solution adapted to each patient. It comprises a surgeon-approved 3D-printed patient-specific surgical guide and patient-specific plate with a small titanium wedge to confirm the 3D angular correction.

[Bodycad BC Fine Osteotomy](#) includes calibrated drill guides with drill stop technology to protect posterior neurovascular structures, multi-planar 3D planning to increase accuracy and is delivered as a procedure in a box to drive efficiency in the O.R.

Contact:

Laurence Robichaud
2035 rue du Haut-Bord
Québec, QC, Canada
G1N 4R7
(418) 928-3391
lrobichaud@bodycad.com